



CAREASSIST® BED

ADVANTA™ 2 BED

VERSACARE® BED

**Feature**

<b>IntelliDrive® Powered Transport:</b> Allows a single caregiver to easily transport patient in bed	Optional	Optional	Optional
<b>Advanced Three-Mode Bed Exit System (Patient Position Monitor):</b> Bed exit system with three different sensitivities	Optional	Optional	Optional
<b>Emergency CPR:</b> Release lever lowers the head section with dampened descent	+	+	+
<b>Brake Alarms:</b> Alerts caregiver when the bed is plugged in and the brake is not set	+	+	+
<b>Chair Position:</b> Caregivers can place patient in a chair position with the touch of one button; bed remains at a low height	+	+	+
<b>Safety Alerts at Foot End of Bed:</b> Provides measures to monitor protocols for bed status (brakes not set, bed exit, bed not in lowest position, siderails down); lights flash an amber color alert when bed is not set in any of the programmed protocols	Optional	Optional	Optional
<b>Integrated Scale</b>	Optional	Optional	Optional
<b>Battery Backup:</b> Backup for all bed articulation (high/low position, head-of-bed, knee section)	+	+	+
<b>Point-of-Care® Siderail Controls:</b> All caregiver controls are conveniently located in a central location, on both sides of the bed, where patient and caregiver interaction occurs	+	+	+
<b>Line-of-Site® Feature:</b> Head-of-bed and Trendelenburg indicators integrated into siderails	+	+	+
<b>HandsFree® Controls:</b> Foot pedals activate high/low and head up/down so your hands are available for your patient			+
<b>Electric Vascular Foot:</b> Eliminates the need for the nursing staff to manually lift the patient's legs during a cardiovascular event	+	+	
<b>Drop Down Siderails:</b> One-handed single action siderail operation	+	+	+
<b>Boost® One-Button Function:</b> Puts the bed in the ideal position for patient repositioning			Optional
<b>Built-in Bed Extender:</b> Easily adds length to the bed	+	+	+
<b>FlexAfoot™ Mechanism:</b> Allows caregiver to customize via siderail controls the overall length of the bed and mattress, in any bed position (including chair position) to fit the patient			+
<b>Ability to send data to Nurse Call, EMR or other Clinical Applications</b>	Optional	Optional	Optional
<b>Fully Integrated Surface:</b> Surface automatically adjusts to the patient's weight, size, movement and to the head-of-bed articulation, to improve patient comfort			Optional
<b>Advanced Microclimate® Technology:</b> Helps reduce heat and moisture at the patient sleep surface interface			Optional
<b>Status Board Nurse Call Display via smart beds</b>	Optional	Optional	Optional



**Feature**

**TOTALCARE® BED**

**TOTALCARE®  
BARIATRIC PLUS BED**

<b>IntelliDrive® Powered Transport:</b> Allows a single caregiver to easily transport patient in bed	Optional	Optional
<b>Advanced Three Mode Bed Exit System (Patient Position Monitor):</b> Bed exit system with three different sensitivities	Optional	Optional
<b>Emergency CPR:</b> One-step CPR flattens the frame and max inflates the air surface	+	+
<b>Brake Alarms:</b> Alerts caregiver when the bed is plugged in and the brake is not set	+	+
<b>Seat Deflate:</b> Seat and foot deflate offers low bed height for easier exit	Optional	Optional
<b>Turn Assist:</b> Uses the bed's pressure redistribution surface to help turn the patient, making it easier to perform linen changes and wound checks	Optional	Optional
<b>CLRT and Percussion &amp; Vibration:</b> Integrated continuous lateral rotation therapy and percussion and vibration modules	Optional	Optional
<b>Preliminary Tilt Table:</b> Extends Reverse Trendelenburg to achieve 20° to gradually acclimate your patient to the upright position	+	+
<b>FullChair® Position:</b> One-button FullChair® position means patients are up in upright FullChair® position with dependent leg positioning for "up in chair" orders	+	+
<b>Chair Egress with Stand Assist:</b> Use the bed to help stand heavier or less sturdy patients	+	+
<b>Integrated Scale</b>	Optional	Optional
<b>Battery Backup:</b> Backup for all bed articulation (high/low position, head-of-bed, knee section)	+	+
	<i>(Dual Electric and Manual Foot Pedal)</i>	<i>(Dual Electric and Manual Foot Pedal)</i>
<b>Point-of-Care® Siderail Controls:</b> All caregiver controls are conveniently located in a central location, on both sides of the bed, where patient and caregiver interaction occurs	+	+
<b>Head-of-Bed Compliance:</b> Multiple digital head-of-bed angle indicators and a head-of-bed alarm make it easier to maintain head-of-bed compliance	+	+
	<i>(Able to store up to 8 days of statistics)</i>	<i>(Able to store up to 8 days of statistics)</i>
<b>Drop Down Siderails:</b> One handed single action siderail operation	+	+
<b>Boost® One-Button Function:</b> Puts the bed in the ideal position for patient repositioning	+	+
<b>FlexAfoot™ Mechanism:</b> Allows caregiver to customize via siderail controls the overall length of the bed and mattress, in any bed position (including chair position) to fit the patient	+	+
<b>Ability to send data to Nurse Call, EMR or other Clinical Applications</b>	Optional	Optional
<b>Fully Integrated Surface:</b> Surface automatically adjusts to the patient's weight, size, movement and to the head-of-bed articulation, to improve patient comfort	Optional	Optional <i>(40" Standard Width)</i>
<b>Advanced Microclimate® technology:</b> Helps reduce heat and moisture at the patient sleep surface interface	Optional	Optional
<b>Status Board Nurse Call Display via smart beds</b>	Optional	Optional

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2014 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

ORDER NUMBER 173672 rev 2 04-MAR-2014 ENG – US

Zinc Number US-AC-1213-0662

[www.hill-rom.com](http://www.hill-rom.com)

USA 800-445-3730  
CANADA 800-267-2337

Enhancing outcomes for patients and their caregivers:

**Hill-Rom**